

Table of Contents

Cell Line / Primary Cell Type	Page
1. K562 (Tier 1)	2
2. GM12878 (Tier 1)	4
3. HeLa S3 (Tier 2)	6
4. HepG2 (Tier 2)	8
5. HUVEC – Human Umbilical Vein Cells (Tier 2)	10
6. NHEK – Normal Human Epidermal Keratinocytes (Tier 2)	13
7. GM06990 (Tier 3)	16
8. SK-N-SH (Tier 3)	18
9. Ovcara-3 (Tier 3)	20
10. BJ-tert (Tier 3)	22
11. Caco-2 (Tier 3)	24
12. HEK293 (Tier 3)	26
13. MCF-7 (Tier 3)	28
14. NTERA-2 (Tier 3)	30
15. Primary Th1 T cells	32

SOP: Propagation of K562
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

K562 may be ordered from ATCC as a frozen ampoule.

Name: K562-chronic myelogenous leukemia
ATCC #: CCL-243

Notes:

This cell line grows in suspension.

Materials List

1. RPMI 1640 with 2mM L-glutamine (cellgro Cat# 10-040-CM)
2. Fetal Bovine Serum (cellgro Cat# 35-016-CV)
3. Sodium Pyruvate (cellgro Cat# 25-000-CI)
4. T225 & T25 culture flasks
5. Graduated pipets (1, 5, 25mL)
6. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
7. Hemocytometer
8. Micropipet w/ P20 tips
9. Microscope

Growth Media for K562

RPMI 1640 with 2mM L-glutamine
Sodium Pyruvate 10mM
10% FBS
Pen-Strep (1x)

Procedure

A. Receipt of Frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath.
- 3) Transfer thawed cells to a T25 flask with 10ml of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) The take cell count and spin down cells, 500g for 5 minutes, then decant old media
- 6) Re-suspend cells in warm fresh media at a volume to yield a density of 2x10⁵ cells/ml.

B. Sub-culture

- 1) Take cell counts every 48 to 72 hours.
- 2) Maintain cell density between 1x10⁵ and 1x10⁶ cells/ml.
- 3) Add fresh warm media when appropriate.
- 4) Record each subculture event as a passage.
- 5) Cells can be spun down, 500g for 5 minutes, rinsed with 1x PBS and re-suspend in a smaller volume of warm growth media when appropriate

C. Maintenance

- 1) Change media as cell density requires.

D. Harvest

- 1) Pass cells until desired cell number is achieved
- 2) Spin down cells and rinse with 1x PBS as described in `sub-culturing`

SOP: Propagation of GM12878
Date modified: 10/29/2008
Modified by: Jeff Goldy/M. Dorschner

Ordering Information

GM12878 may be ordered from Coriell Cell Repositories. Proliferating cells are shipped in a T25 flask with 10-20ml of media.

To order starter cultures:

Name/Catalogue #: GM12878

Notes:

This cell line grows in suspension and should be maintained at a density between 2×10^5 cells/ml and 1×10^6 cells/ml.

Materials List

1. RPMI 1640 with 2mM L-glutamine (cellgro Cat# 10-040-CM)
2. Fetal Bovine Serum (cellgro Cat# 35-016-CV)
3. T225 culture flasks
4. Graduated pipets (1, 5, 25, 50mL)
5. Penicillin-Streptomycin Solution, 100X (Cellgro, Cat#300-002CI)
6. Hemocytometer
7. Micropipet w/ P20 tips
8. Microscope
9. Freezing medium (growth medium containing 6% DMSO)

Growth Medium for GM12878

RPMI 1640 with 2mM L-glutamine
15% FBS
Pen-Strep (1X)

Procedure

A. Receipt of proliferating cells and generation of seed stocks

- 1) Equilibrate unopened T25 flask overnight in 37°C, 5% CO₂ humidified incubator to allow cells to recover.
- 2) Cells should be counted the next day and split to achieve a cell density of 200,000-500,000 cells/ml.
- 3) Cells should be incubated in upright flasks with vented or loose caps.
- 4) Upon reaching the desired number, cells should be spun down, rinsed with 1X PBS, resuspended in freezing medium.
- 5) Cells are dispensed into cryovials (2 million per aliquot) and frozen in a -80°C isopropanol bath overnight.
- 6) Cryovials are transferred the next day to liquid nitrogen for long term storage.

B. Sub-culture and Maintenance

- 1) Maintain culture at a cell density between 2×10^5 and 1×10^6 cells/ml.
- 2) Cells will either need to be fed again after 3-4 days or split depending on the cell density. Splitting can be performed by centrifuging cells at 500g for 5 minutes, decanting growth medium and rinsing in sterile 1X PBS. Cells should then be resuspended in fresh growth medium to achieve a density 2×10^5 and 1×10^6 cells/ml.

C. Harvest

- 1) Pass cells until the desired number of cells is reached.
- 2) Spin down and rinse cells as described above in Sub-culture and maintenance.

SOP: Propagation of HeLa S3
Date modified: 10/29/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

HeLa S3 may be ordered from ATCC as a frozen ampoule.

Name: HeLa S3, cervical carcinoma
ATCC #: CCL-2.2

Notes:

This is an adherent cell line. We use DMEM in place of ATCC recommended F-12K medium.

Materials List

1. DMEM with 2mM L-glutamine (cellgro Cat# 10-013-CM)
2. Fetal Bovine Serum (cellgro Cat# 35-016-CV)
3. T75 & T225 culture flasks
4. Graduated pipets (1, 5, 25mL)
5. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
6. Hemocytometer
7. Micropipet w/ P20 tips
8. Microscope

Growth Media for HeLa S3

DMEM with 2mM L-glutamine
10% FBS
1x Pen-Strep

Procedure

A. Receipt of Frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37C water bath
- 3) Transfer thawed cells to a T75 flask with 40ml of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 3 minutes (4°C)
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:4 to 1:8 cell split as needed
- 9) Record each subculture event as a passage

C. Maintenance

- 1) Change media the day after seeding and every 72 hours thereafter.
Use 50ml of media per T225

D. Harvest

- 1) Do not use cells that have been passed more than 25 times
- 2) Remove cells from flasks according to protocol described above under 'sub-culturing'
- 3) Examine viability using trypan blue staining.

SOP: Propagation of HepG2
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

HepG2 can be ordered from ATCC as a frozen ampoule.

Name: HepG2, hepatocellular carcinoma
ATCC #: HB-8065

Notes:

This is an adherent cell line.

Materials List

1. MEM with 2mM L-glutamine (Cellgro Cat# 10-010-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Sodium pyruvate (Cellgro Cat# 25-000-CI)
5. Non-essential amino acids (Cellgro Cat# 25-025-CI)
6. T75 & T225 culture flasks
7. Graduated pipets (1, 5, 25mL)
8. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
9. Hemocytometer
10. Micropipet w/ P20 tips
11. Microscope

Growth Medium for HepG2

MEM with 2mM L-glutamine
Non-essential amino acids
Sodium Pyruvate 1mM
10% FBS
Pen-Strep (1X)
Sodium Bicarbonate 1.5g/L

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37 °C water bath
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C).
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:4 to 1:6 cell split as needed.
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and every 2-3 days thereafter.
Use ~50mLs of medium per T225 flasks.

D. Harvest

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Examine viability using trypan blue staining (SOP).

SOP: Propagation of primary Human Umbilical Vein Endothelial Cells
(HUVECs; Lonza Biosciences)
Date modified: 8/27/08
Modified by: M. Dorschner (UW)
L. Dillon (NHGRI)

Ordering Information

Primary Human Umbilical Vein Endothelial Cells (HUVECs) may be ordered either as frozen ampules or as starter cultures. The former contain $\sim 0.5-1 \times 10^5$ cells; the latter are initiated at Lonza and sent in a T225 flask containing $6-7 \times 10^6$ cells.

For all orders, provide (1) Reservation #; (2) Contract/quotation #; (3) Individual (Lot #); and (4) Item #s, as follows:

Reservation number: ~~RZ-495718~~ 3122124 (Updated 7/22/08)
Contract number: P101416
Individual H1: Lot #7F3239 - Male, Caucasian,
154 amps available, 0 amps available (7/22/08)

Individual H2: Lot #7F3771 - Male, African American
144 amps available, 38 amps available (7/22/08)

To order frozen ampules + media:

Name: HUVEC – Umbilical Vein Endo Cells
Item #: CC-2517 (HUVEC in EGM® - Cryopreserved ampule)
CC-3162 (EGM-2 BulletKit = CC-3156 + CC-4176)

To order starter cultures:

Name: HUVEC – Umbilical Vein Endo Cells
Item #: CC2501T225 (HUVECs in EGM® T-225 Flask)
CC-3162 (EGM-2 BulletKit = CC-3156 + CC-4176)

Notes:

The number of BulletKits purchased depends on the target number of cells to be generated. A rule of thumb is 10 BulletKits for every initial T225 flask of cells. It is strongly recommended to purchase all of the media that will be required for a complete expansion series (see below), since media supply may be erratic.

Materials List

1. Cell-type specific medium (BulletKits – Lonza Biosciences)
2. T225 culture flasks
3. Graduated pipets (1, 5, 25mL)
4. Pen-strep solution (if required; Lonza typically supplies antibiotics)
5. Hemocytometer
6. Micropipet w/ P20 tips
7. Microscope

Procedure

A. Receipt of proliferating cells

- 1) Equilibrate for 3-4 hours in 37°C, 5% CO₂ humidified incubator.
- 2) Remove shipping medium. Replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 3 minutes (4°C)
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Count cells with hemocytometer.
- 9) Add warmed medium to flasks.
- 10) Seed flasks at **5,000 cells/cm²**
- 11) Record each subculture event as a passage

C. Maintenance


- 1) Change media the day after seeding and every OTHER day thereafter.
- 2) Increase media volume as confluency increases (volumes assume the use of
- 3) T225 flasks):
 - a. 25 % = 1mL/5 cm²
 - b. 25-45% = 1.5mL/ 5 cm²
 - c. 45%+ = 2mL/ 5 cm²
- 4) Per the above an exemplary schedule might be:
 - a. day 1, plate into T225: use 50 mls of media.
 - b. day 2, change media, use 50 mls of media
 - c. day 4, change media, use 100 mls of media (if confluency is >50%)
 - d. day 6, change media, use 100 mls of media (or harvest if ready)
 - e. day 7 or 8 (harvest when cells reach 6 x 10⁶ cells/flask)

D. Harvest

- 1) Pass cells 3-4 times until the desired cell number is achieved (primary cells will senesce after 4-5 passages).
- 2) Remove cells from flasks according to protocol described above under 'subculturing'
- 3) Examine viability using trypan blue staining (SOP)

Exemplary Expansion

The diagram below illustrates an exemplary expansion of HUVECs from a Lonza starter culture:



- The initial T225 flask received from Lonza will have $\sim 6 \times 10^6$ cells; this will then be split and seeded at $\sim 3,500$ cells/cm²; each new T225 flask will therefore be seeded with $\sim 750K$ cells.
- The initial flask will yield up to 7-8 daughter flasks depending on how large of an expansion is targeted.
- Once these flasks have reached the target density again, they can be split and seeded into up to 40 flasks.
- The 40 granddaughter flasks will each yield $\sim 6 \times 10^6$ cells, providing a total theoretical yield of 2.5×10^8 cells.

Media requirements: Each flask will require ~ 50 mL of medium with additional medium for feedings during the doubling/expansion process.

SOP: Propagation of primary keratinocytes (NHEK; Lonza Biosciences)
Date modified: 7/22/08
Modified by: M. Dorschner (UW)
L. Dillon (NHGRI)

Ordering Information

Normal Human Epidermal Keratinocytes (NHEKs) may be ordered either as frozen ampules or as starter cultures. The former contain $\sim 0.5-1 \times 10^5$ cells; the latter are initiated at Lonza and sent in a T225 flask containing $6-7 \times 10^6$ cells.

For all orders, provide (1) Reservation #; (2) Contract/quotation #; (3) Individual (Lot #); and (4) Item #s, as follows:

Reservation number: ~~RZ-495718~~ 3122124 (Updated 7/22/08)
Contract number: P101416
Individual K1: Lot #4F1155J - Female, African American
~~178 amps available~~, 51 amps available (7/22/08)

Individual K2: Lot #7F4307 - Female, Caucasian
~~179 amps available~~; 139 amps available (7/22/08)

To order frozen ampules + media:

Name: NHEK – Adult Keratinocytes
Item #: CC-2501 (NHEK in KGM® - Cryopreserved ampule)

CC-3111 (KGM BulletKit = CC-3101 + CC-4131)

To order starter cultures:

Name: NHEK – Adult Keratinocytes
Item #: CC2501T225 (NHEKs in KGM® T225 Flask)
CC-3111 (KGM BulletKit = CC-3101 + CC-4131)

Notes:

The number of BulletKits purchased depends on the target number of cells to be generated. A rule of thumb is 10 BulletKits for every initial T225 flask of cells. It is strongly recommended to purchase all of the media that will be required for a complete expansion series (see below), since media supply may be erratic.

Materials List

1. Cell-type specific medium (BulletKits – Lonza Biosciences)
2. T225 culture flasks
3. Graduated pipets (1, 5, 25mL)
4. Pen-strep solution (if required; Lonza typically supplies antibiotics)
5. Hemocytometer
6. Micropipet w/ P20 tips
7. Microscope

Procedure

A. Receipt of proliferating cells

- 1) Equilibrate for 3-4 hours in 37°C, 5% CO₂ humidified incubator.
- 2) Remove shipping medium. Replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 3 minutes (4°C)
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Count cells with hemocytometer.
- 9) Add warmed medium to flasks.
- 10) Seed flasks at **3,500 cells/cm²**
- 11) Record each subculture event as a passage

C. Maintenance


- 1) Change media the day after seeding and every OTHER day thereafter.
- 2) Increase media volume as confluency increases (volumes assume the use of
- 3) T225 flasks):
 - a. 25 % = 1mL/5 cm²
 - b. 25-45% = 1.5mL/ 5 cm²
 - c. 45%+ = 2mL/ 5 cm².
- 4) Per the above an exemplary schedule might be:
 - a. day 1, plate into T225: use 50 mls of media.
 - b. day 2, change media, use 50 mls of media
 - c. day 4, change media, use 100 mls of media (if confluency is >50%)
 - d. day 6, change media, use 100 mls of media (or harvest if ready)
 - e. day 7 or 8 (harvest when cells reach 6 x 10⁶ cells/flask

D. Harvest

- 1) Pass cells 3-4 times until the desired cell number is achieved (primary cells will senesce after 4-5 passages).
- 2) Remove cells from flasks according to protocol described above under 'subculturing'
- 3) Examine viability using trypan blue staining (SOP TP-7)

Exemplary Expansion

The diagram below illustrates an exemplary expansion of NHEKs from a Lonza starter culture:



- The initial T225 flask received from Lonza will have ~6 x 10⁶ cells; this will then be split and seeded at ~3,500 cells/cm²; each new T225 flask will therefore be seeded with ~750K cells.
- The initial flask will yield up to 7-8 daughter flasks depending on how large of an expansion is targeted.
- Once these flasks have reached the target density again, they can be split and seeded into up to 40 flasks.
- The 40 granddaughter flasks will each yield ~6 x 10⁶ cells, providing a total theoretical yield of 2.5x10⁸ cells.

Media requirements: Each flask will require ~50mL of medium with additional medium for feedings during the doubling/expansion process.

SOP: Propagation of GM06990
Date modified: 10/29/2008
Modified by: Jeff Goldy/M. Dorschner

Ordering Information

GM06990 may be ordered from Coriell Cell Repositories. Proliferating cells are shipped in a T25 flask with 10-20ml of media.

To order starter cultures:

Name/Catalogue #: GM06990

Notes:

This cell line grows in suspension and should be maintained at a density between 2×10^5 cells/ml and 1×10^6 cells/ml.

Materials List

1. RPMI 1640 with 2mM L-glutamine (cellgro Cat# 10-040-CM)
2. Fetal Bovine Serum (cellgro Cat# 35-016-CV)
3. T225 culture flasks
4. Graduated pipets (1, 5, 25, 50mL)
5. Penicillin-Streptomycin Solution, 100X (Cellgro, Cat#300-002CI)
6. Hemocytometer
7. Micropipet w/ P20 tips
8. Microscope
9. Freezing medium (growth medium containing 6% DMSO)

Growth Medium for GM06990

RPMI 1640 with 2mM L-glutamine

15% FBS

Pen-Strep (1X)

Procedure

A. Receipt of proliferating cells and generation of seed stocks

- 1) Equilibrate unopened T25 flask overnight in 37°C, 5% CO₂ humidified incubator to allow cells to recover.
- 2) Cells should be counted the next day and split to achieve a cell density of 200,000-500,000 cells/ml.
- 3) Cells should be incubated in upright flasks with vented or loose caps.
- 4) Upon reaching the desired number, cells should be spun down, rinsed with 1X PBS, resuspended in freezing medium.
- 5) Cells are dispensed into cryovials (2 million per aliquot) and frozen in a -80°C isopropanol bath overnight.
- 6) Cryovials are transferred the next day to liquid nitrogen for long term storage.

B. Sub-culture and Maintenance

- 1) Maintain culture at a cell density between 2×10^5 and 1×10^6 cells/ml.
- 2) Cells will either need to be fed again after 3-4 days or split depending on the cell density. Splitting can be performed by centrifuging cells at 500g for 5 minutes, decanting growth medium and rinsing in sterile 1X PBS. Cells should then be resuspended in fresh growth medium to achieve a density 2×10^5 and 1×10^6 cells/ml.

C. Harvest

- 1) Pass cells until the desired number of cells is reached.
- 2) Spin down and rinse cells as described above in Sub-culture and maintenance.

SOP: Propagation of Sk-N-SH
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

SK-N-SH can be ordered from ATCC as a frozen ampule.

Name: Sk-N-SH, Neuroblastoma
ATCC #: HTB-11

Notes:

This is an adherent cell line. Cells are differentiated with retinoic acid (ATRA) for 48 hours prior to harvesting.

Materials List

1. RPMI 1640 with 2mM L-glutamine (Cellgro Cat# 10-040-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Pyruvate (Cellgro Cat# 25-000-CI)
4. T225 culture flasks
5. *all trans*-Retinoic acid (Sigma, Cat #R2625)
6. Graduated pipets (1, 5, 25mL)
7. Penicillin-Streptomycin Solution, 100X (Cellgro, Cat#300-002CI)
8. Hemocytometer
9. Micropipet w/ P20 tips
10. Microscope

Growth Medium for SK-N-Sh

RPMI 1640 with 2mM L-glutamine
Sodium Pyruvate 10mM
10% FBS
Pen-strep (1X)

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath
- 3) Transfer thawed cells to a T75 flask with 40ml of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 3 minutes (4°C).
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:3 to 1:8 cell split as needed.
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and 1-2 times per week thereafter.
Use ~50mLs of medium per T225 flask.

D. Differentiation

- 1) Upon reaching the desired cell number, add growth medium containing 6uM *all trans*-retinoic acid. Cells should be cultured for 48 hours in differentiation medium.

E. Harvest

- 1) Do not use cells that have been passed more than 8 times
- 2) Remove cells from flasks according to protocol described above under 'subculturing'
- 3) Examine viability using trypan blue staining (SOP)

SOP: Propagation of Ovc3
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

Ovc3 can be ordered from ATCC as a frozen ampoule.

Name: Ovc3, ovarian adenocarcinoma
ATCC #: HTB-161

Notes:

This is an adherent cell line.

Materials List

1. RPMI-1640 with 2mM L-glutamine (Cellgro Cat# 10-040-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Sodium Pyruvate (Cellgro Cat# 25-000-CI)
5. HEPES buffer (Cellgro Cat# 25-060-CI)
6. T75 & T225 culture flasks
7. Graduated pipets (1, 5, 25mL)
8. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
9. Hemocytometer
10. Micropipet w/ P20 tips
11. Microscope

Growth Medium for Ovc3

RPMI-1640 with 2mM L-glutamine
HEPES Buffer 10mM
Sodium Pyruvate 1mM
20% FBS
Pen-Strep (1X)
Sodium Bicarbonate 1.5g/L
Bovine Insulin 0.01mg/ml

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath.
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C)
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:2 to 1:4 cell split as needed.
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and every 2-3 days thereafter.
Use ~50ml of medium per T225 flask.

D. Harvest

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Let cells grow 48 hours past confluence.
- 4) Examine viability using trypan blue staining (SOP).

SOP: Propagation of BJ-tert
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

BJ-tert can be ordered from ATCC as a frozen ampoule.

Name: BJ-tert, skin fibroblast
ATCC #: CRL-2522

Notes:

This is an adherent cell line.

Materials List

1. MEM with 2mM L-glutamine (Cellgro Cat# 10-010-CM).
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Sodium pyruvate (Cellgro Cat# 25-000-CI)
5. Non-essential amino acids (Cellgro Cat# 25-025-CI)
6. T75 & T225 culture flasks
7. Graduated pipets (1, 5, 25mL)
8. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
9. Hemocytometer
10. Micropipet w/ P20 tips
11. Microscope

Growth Medium for BJ-tert

MEM with 2mM L-glutamine
Non-essential amino acids
Sodium Pyruvate 1mM
10% FBS
Pen-Strep (1X)
Sodium Bicarbonate 1.5g/L

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C).
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:2 to 1:9 cell split as needed.
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and every 2-3 days thereafter.
Use ~50ml of medium per T225 flasks.

D. Harvest

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Examine viability using trypan blue staining (SOP).

SOP: Propagation of Caco-2
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

Caco-2 can be ordered from ATCC as a frozen ampoule.

Name: Caco-2, colorectal adenocarcinoma
ATCC #: HTB-37

Notes:

This is an adherent cell line.

Materials List

1. MEM with 2mM L-glutamine (Cellgro Cat# 10-010-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Sodium Pyruvate (Cellgro Cat# 25-000-CI)
5. Non-essential amino acids (Cellgro Cat# 25-025-CI)
6. T75 & T225 culture flasks
7. Graduated pipets (1, 5, 25mL)
8. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
9. Hemocytometer
10. Micropipet w/ P20 tips
11. Microscope

Growth Medium for Caco-2

MEM with 2mM L-glutamine
Non-essential amino acids
Sodium Pyruvate 1mM
20% FBS
Pen-Strep (1X)
Sodium Bicarbonate 1.5g/L

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath.
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches ~80% confluence ($8 \times 10^4 - 10^5$ cells/cm²).
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C)
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:4 to 1:8 cell split as needed. Seed at a density of 10^4 cells/cm².
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and 1- 2 times per week.
Use ~50mLs of medium per T225 flask.

D. Differentiation

- 1) Upon reaching the desired cell number, cells are grown to confluence. Cells are not harvested until 2 days after confluence to ensure complete differentiation.

E. Harvest

- 2) Do not use cells that have been passed more than 8 times.
- 3) Remove cells from flasks according to protocol described above under 'subculturing'.
- 4) Let cells grow 48 hours past confluence.
- 5) Examine viability using trypan blue staining (SOP).

SOP: Propagation of HEK293
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

Hek293 can be ordered from ATCC as a frozen ampoule.

Name: HEK293, embryonic kidney
ATCC #: CRL-1573

Notes:

This is an adherent cell line.

Materials List

1. MEM with 2mM L-glutamine (Cellgro Cat# 10-010-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Sodium pyruvate (Cellgro Cat# 25-000-CI)
5. Non-essential amino acids (Cellgro Cat# 25-025-CI)
6. T225 culture flasks
7. Graduated pipets (1, 5, 25mL)
8. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
9. Hemocytometer
10. Micropipet w/ P20 tips
11. Microscope

Growth Medium for HEK293

MEM with 2mM L-glutamine
Non-essential amino acids
Sodium Pyruvate 1mM
10% FBS
Pen-strep (1X)
Sodium Bicarbonate 1.5g/L

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath.
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C)
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:2 to 1:4 cell split as needed.
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and every 2-3 days thereafter.
Use ~50mLs of medium per T225 flask.

D. Harvest

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Examine viability using trypan blue staining (SOP).

SOP: Propagation of MCF-7
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

MCF-7 can be ordered from ATCC as a frozen ampoule.

Name: MCF-7, mammary gland, adenocarcinoma
ATCC #: HTB-22

Notes:

This is an adherent cell line.

Materials List

1. MEM with 2mM L-glutamine (Cellgro Cat# 10-010-CM)
2. Fetal Bovine Serum (cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (cellgro Cat# 25-035-CI)
4. Non-essential amino acids (cellgro Cat# 25-025-CI)
5. T75 & T225 culture flasks
6. Graduated pipets (1, 5, 25mL)
7. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
8. Hemocytometer
9. Micropipet w/ P20 tips
10. Microscope

Growth Medium for MCF-7

MEM with 2mM L-glutamine
Non-essential amino acids
10% FBS
Pen-Strep (1X)
Sodium Bicarbonate 1.5g/L

Procedure

A. Receipt of frozen cells and starting cell cultures

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampule in 37°C water bath.
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Decant medium.
- 3) Wash cells with warm 1X PBS.
- 4) Add 8mLs of Accutase and return to incubator for 10-15 minutes.
- 5) Immediately remove cells and pellet at 500 xg for 5 minutes (4°C).
- 6) Wash cells 2X with 1X PBS.
- 7) Gently re-suspend cell pellet in warm medium.
- 8) Perform 1:3 to 1:6 cell split as needed.
- 9) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and every 2-3 days thereafter.
Use ~50mLs of medium per T225 flask.

D. Harvest

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Examine viability using trypan blue staining (SOP).

SOP: Propagation of NTERA-2
Date modified: 9/5/2008
Modified by: J. Goldy/M. Dorschner

Ordering Information

NTERA-2 can be ordered from ATCC as a frozen ampoule.

Name: NTERA-2, testis-malignant pluripotent embryonal carcinoma
ATCC #: CRL-1973

Notes:

This is an adherent cell line. Cells must be scraped for subculturing and harvesting.

Materials List

1. DMEM with 2mM L-glutamine (Cellgro Cat# 10-013-CM)
2. Fetal Bovine Serum (Cellgro Cat# 35-016-CV)
3. Sodium Bicarbonate (Cellgro Cat# 25-035-CI)
4. Cell Scraper (Falcon Cat# 353087)
5. T75 & T225 culture flasks
6. Graduated pipets (1, 5, 25mL)
7. Penicillin-Streptomycin Solution (100X) (Cellgro Cat# 30-002-CI)
8. Hemocytometer
9. Micropipet w/ P20 tips
10. Microscope

Growth Medium for NTERA-2

DMEM with 2mM L-glutamine
10% FBS
Pen-strep (1X)
Sodium Bicarbonate 1.5g/L

Procedure

A. Receipt of frozen cells and starting cell cultures.

- 1) Immediately place frozen cells in liquid nitrogen storage incubator.
- 2) Quickly thaw ampoule in 37°C water bath.
- 3) Transfer thawed cells to a T75 flask with 40mLs of warm growth media.
- 4) Allow cells to recover over night in 37°C, 5% CO₂ humidified incubator.
- 5) Pour off medium the next day, replace with fresh medium and return to incubator.

B. Sub-culture

- 1) Propagate cells until density reaches 70-80% confluence.
- 2) Scrape cells from each flask using cell scraper.
- 3) Immediately remove cells each flask and pellet at 500 xg for 5 minutes (4°C).
- 4) Gently re-suspend cell pellet in warm medium.
- 5) Perform 1:4 to 1:8 cell split as needed. New subcultures should be seeded at $\geq 1.5 \times 10^7$ viable cells per 225 cm² flask.
- 6) Record each subculture event as a passage.

C. Maintenance

- 1) Change media the day after seeding and every 2-3 days thereafter. Cultures should be maintained at high density. Use ~50mLs of medium per T225 flask.

D. Harvest

- 1) Do not use cells that have been passed more than 8 times.
- 2) Remove cells from flasks according to protocol described above under 'subculturing'.
- 3) Examine viability using trypan blue staining (SOP).

SOP: Propagation of Th1
Date modified: 11/17/2008
Modified by: M. Dorschner

Source Information:

Cells are procured from primary pheresis of a single normal subject.

Notes:

Th1 T-cell subset is purified and expanded in primary culture.

Materials List

1. Naïve CD4+ T cell isolation kit (Miltenyi Biotech, Cat # 130-091-894)
2. autoMACS Separator (Miltenyi Biotech)
3. Aim V medium (Invitrogen, Cat # 087-0112DK)
4. AB serum (Lonza Bioscience, Cat # 14-490E)
5. Anti-CD3/Anti-CD28 coated beads (Dyna/Invitrogen, Cat #111-31D)
6. Human IL-2 (R&D Systems)
7. Human IL-12 (R&D Systems)
8. Anti-Human IL-4 (eBioscience)
9. Bio-Plex Human Cytokine Th1/Th2 Panel (Bio-Rad, Cat # 171-A11081)
10. T25 & T225 culture flasks
11. Graduated pipets
12. Hemocytometer
13. Phorbol 12-myristate 13-acetate (Sigma, Cat # P1585)
14. Ionomycin (Sigma, Cat # I3909)

Th1 Polarization Medium

Aim V medium
2% AB serum
Human IL-2 (50 IU/mL)
Human IL-12 (10 ng/mL)
Anti-Human IL-4 (5 ug/mL)

Procedure

A. Isolation of naïve CD4+ T cells

- 1) Isolate naïve CD4+ T cells by negative selection using the Naïve CD4+ T Cell Isolation Kit according to manufacturer's recommendations.

B. Stimulation and polarization of cells

- 1) Resuspend naïve CD4⁺ T cells (93% \geq CD4⁺CD45RA⁺) in Aim V medium containing 2% serum.
- 2) Stimulate cells with anti-CD3 and anti-CD28 coated beads in polarizing medium.
- 3) Expand cells in culture for 7-10 days.
- 4) If needed, stimulate cells with PMA (2.5 ug/mL) and Ionomycin (500 uM) to divide.

C. Confirmation of Th1 polarization

- 1) Assay supernatant with Bio-Plex Human Cytokine Th1/Th2 panel according to manufacturer's protocol.

D. Harvest

- 1) Pellet cells by centrifugation and wash cells in 1X PBS.
- 2) Examine viability using trypan blue staining.